

# PRODUCT NOTE

## FMS-CARBON DIOXIDE Headspace Analyzer



The FMS-Carbon Dioxide is a non-destructive headspace carbon dioxide analyzer from LIGHTHOUSE, the industry leader in laser-based headspace analysis. The compact benchtop instrument utilizes a patented laser absorption technique developed with funding from the Food and Drug Administration. This rapid and versatile technology addresses a wide range of applications that span the full product life cycle.

## APPLICATIONS

- Container closure integrity testing of frozen product stored on dry ice for transport
- General container closure integrity testing when the test method uses carbon dioxide as a tracer gas
- IPC monitoring of carbon dioxide levels during the filling of product purged with carbon dioxide in the headspace
- Microbial growth detection in media vials

## KEY FEATURES

- Non-destructive, quantitative measurement method
- High-sensitivity signal analysis delivers an accurate measurement in seconds
- Custom change parts provide consistent positioning of sample for optimal measurement across a wide range of container types and sizes
- NIST-certified carbon dioxide standards for calibration and verification to ensure accurate results
- Easy-to-use hardware and software requires minimal user training
- Full validation package and 21 CFR Part 11 compliant software

## SYSTEM SPECIFICATIONS

- Measurement range: 0.0 – 1.0 atm
- Measurement time: 0.5 – 5.0 seconds
- Sample type: syringe, ampoule, vial, bottle
- Sample size: 6.0 – 86.0 mm in diameter (1ml syringe to 200ml bottle)
- Dimensions: 30.5 x 30.5 x 29.2 cm
- Weight: 13.6 kg
- Power requirements: 110 - 240 VAC, 50/60 Hz, 60W
- Interface: PC
- Safety Standards: IEC/EN 61010, 61326, 60825; US CDRH 21 CFR 1040; Declaration of Conformity available

## APPLICATION SPOTLIGHT

The versatility of the headspace carbon dioxide measurement enables the collection of analytical data for a wide range of applications, providing an opportunity for process optimization and improvements to



product quality. The following study results highlight examples of data that can be generated using the FMS-Carbon Dioxide system.

## CONTAINER CLOSURE INTEGRITY TESTING

The use of carbon dioxide as a tracer gas provides a high-sensitivity leak detection method. After using a LIGHTHOUSE CCI Test Vessel to expose vials to carbon dioxide, the positive controls with idealized 5-micron defects were easily identified by the significant increase in headspace carbon dioxide (Figure 1). Initial levels near zero increased to almost 1 atm (760 torr) of carbon dioxide in the leaking positive control vials.




Figure 1: Container Closure Integrity Testing. After exposure to carbon dioxide in the CCI Test Vessel, positive control defects are easily identified by the increase in headspace carbon dioxide.

## CCI TESTING ON DRY ICE

The loss of rubber stopper elasticity under deep cold storage poses an increased risk of seal integrity failure. A leak that develops in dry ice storage will allow the ingress of the surrounding carbon dioxide gas. After storage on dry ice for seven days, 3 of the 20 product samples in this study exhibited an increase in the headspace carbon dioxide (Figure 2).




Figure 2: CCI Testing on Dry Ice. After storage on dry ice, 3 samples reveal elevated headspace carbon dioxide levels due to loss of seal integrity.

## MICROBIAL GROWTH DETECTION IN MEDIA VIALS

Aerobic microbial growth within a media vial will result in the consumption of oxygen and the release of carbon dioxide. A contaminated media vial will exhibit an increase in carbon dioxide during the standard incubation period, providing a rapid and quantitative method for media fill inspection. This result has been observed in a wide range of microorganisms, including the *B. spizizenii* samples shown here (Figure 3), where rising carbon dioxide and falling oxygen were detected within three days of incubation.




Figure 3: Microbial Growth Detection in Media Vials. Headspace oxygen (solid-line) and carbon dioxide (dashed-line) concentrations in media vials injected with *B. spizizenii*. Growth is detected after 24 hrs, with headspace conditions changing drastically in the first 72 hours.

Please contact us or visit our website for more information: [www.lighthouseinstruments.com](http://www.lighthouseinstruments.com)

**Corporate Headquarters**

Lighthouse Instruments, LLC  
2020 Avon Court, Suite #2  
Charlottesville, VA 22902  
Telephone: +1 (434) 293-3081  
E-mail: [info@lighthouseinstruments.com](mailto:info@lighthouseinstruments.com)

**European Support Center**

Lighthouse Instruments B.V.  
Science Park 408  
1098 XH Amsterdam, The Netherlands  
Telephone: +31 20 7051 050  
E-mail: [euinfo@lighthouseinstruments.com](mailto:euinfo@lighthouseinstruments.com)